

Statement of Accuracy Zinc Vacu-vials® Kit, Cat. No. K-9923

This statement of accuracy is based on laboratory tests performed under ideal testing conditions using standards of known concentration prepared in deionized water.

| Instrument Platform | Test Range, ppm | Accuracy (expressed as ppm at concentrations across the test range) | | | |
|---------------------|--------------------|---|---------|---------|----------|
| | | 0 ppm | 1.5 ppm | 3.8 ppm | 11.3 ppm |
| Spectrophotometer | 0 - 15.0 | ≤ 0.4 | +/- 0.5 | +/- 0.8 | +/- 1.1 |
| V-2000 | 0 - 15.0 | ≤ 0.4 | +/- 0.5 | +/- 0.8 | +/- 1.1 |
| V-3000 Series | 0 - 15.0 | ≤ 0.4 | +/- 0.5 | +/- 0.8 | +/- 1.1 |

Updated 12/01/2020